

1 Second Testing

- QUICK
- SIMPLE
- ACCURATE

4 IN 1

- BODY TEMPERATURE
- SURFACE TEMPERATURE
- HIGH TEMPERATURE ALARM
- MEMORIZE 32 TEST DATA

TECHNICAL SPECIFICATIONS

1.Normal using condition:

Ambient temperature:10°C~40°C (50°F~104°F)

Relative humidity:≤85%

2.Batteries: DC 3V (2 pcs AA batteries)

3.Unit size: 155 x 100 x 40 mm (L x W x H)

4.Temperature display resolution: 0.1°C (0.1°F)

5.Measuring range:

In body mode: 32.0°C~43.0°C (89.6°F~109.4°F)

Under body mode, there is three color backlights:

Green color backlit:

≤37.3°C (99.1°F), means normal temperature.

Orange color backlit:

37.4°C~37.9°C (99.3°F~100.2°F), means low fever.

Red color backlit:

≥38°C (100.4°F), means high fever.

In surface temp mode: 0°C~60°C (32°F~140°F)

In room mode: 0°C~40°C(32°F~104°F)

6.Precision:

32.0°C~34.9°C (89.6°F~94.8°F) ±0.3°C(±0.6°F)

35.0°C~42.0°C (95°F~107.6°F) ±0.2°C(±0.4°F)

42.1°C~43.0°C (107.8°F~109.4°F) ±0.3°C(±0.6°F)

7.Consumption: ≤300mW

8.Accuracy: ±0.3°C (0.6°F)

9.Measuring distance: 3cm ~ 5cm (1.2in ~ 2in))

10.Automatic power-off: <30 secs

11.Memory: 32 sets

12. Measuring position: Forehead skin



Please read the product manual carefully or purchase and use under the guidance of medical staff
For details on contraindications or precautions